

For Research Use Only

NOXA2/p67phox Recombinant monoclonal antibody, PBS Only

Catalog Number: 88418-1-PBS



Basic Information

Catalog Number:

88418-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8691

GenBank Accession Number:

BC001606

GeneID (NCBI):

4688

UNIPROT ID:

P19878

Full Name:

neutrophil cytosolic factor 2

Calculated MW:

60 kDa

Observed MW:

67 kDa

Purification Method:

Protein A purification

CloneNo.:

260956B9

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

NOXA2/p67phox is also named as NCF2 (Neutrophil cytosol factor 2) and belongs to the NCF2/NOXA1 family. It is a membrane bound cytochrome b558 required for activation of the latent NADPH oxidase. NCF2 is an essential component of the multi-protein NADPH oxidase enzyme in phagocytic leukocytes, as well as in certain non-phagocytic cells. This protein has 4 isoforms produced by alternative splicing.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88418-1-RR (NCF2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88418-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed THP-1 cells using NOXA2/p67phox antibody (88418-1-RR, Clone: 260956B9) at dilution of 1:500 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 88418-1-PBS in a different storage buffer formulation.