

For Research Use Only

Gamma Actin Recombinant monoclonal antibody, PBS Only

Catalog Number: 88421-1-PBS



Basic Information

Catalog Number:

88421-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1728

GenBank Accession Number:

BC009848

GeneID (NCBI):

71

UNIPROT ID:

P63261

Full Name:

actin, gamma 1

Calculated MW:

42 kDa

Observed MW:

42-45 kDa

Purification Method:

Protein A purification

CloneNo.:

260974G7

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88421-1-RR (ACTG1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88421-1-PBS in a different storage buffer formulation.