

CHAF1A Recombinant monoclonal antibody

Catalog Number: 88451-1-RR

Basic Information

Catalog Number: 88451-1-RR	GenBank Accession Number: BC067093	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 10036	CloneNo.: 260852B11
Source: Rabbit	UNIPROT ID: Q13111	Recommended Dilutions: WB: 1:500-1:2000
Isotype: IgG	Full Name: chromatin assembly factor 1, subunit A (p150)	
Immunogen Catalog Number: AG10536	Calculated MW: 956 aa, 107 kDa	
	Observed MW: 150 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : MOLT-4 cells, K-562 cells
Species Specificity: human	

Background Information

Chromatin assembly factor 1 (CAF1) is the only histone chaperone known to assemble histones H3 and H4 onto newly synthesized DNA both in vitro and in vivo [PMID:17065558]. The 938 amino acid multidomain p150 (CHAF1A) binds via its C-terminal third to p60, which is an essential step for nucleosome assembly because knocking down either subunit disrupts the activity [PMID:14519857]. In addition, CAF1 facilitates DNA synthesis depending on the binding of the N-terminal 31 residues of p150 to the proliferating cell nuclear antigen (PCNA), which acts as a sliding clamp to stimulate the processivity of DNA polymerase [PMID:10648606]. CHAF1A regulates the formation of heterochromatin in mammalian cells during replication and in plants it maintains the transcription of certain subsets of genes. Furthermore, CHAF1A exists in a chromatin-remodeling complex WINAC, which coactivates ligand-induced transactivation function of the vitamin D receptor [PMID:12837248]. CHAF1A protein exists some phosphorylation sites, which may affect its theoretical molecular weight when tested. And a 150 kDa band was recognized (PMID:27445493).

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88451-1-RR (CHAF1A antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.