

For Research Use Only

Anti-Mouse CD86 Rabbit Recombinant Antibody

Catalog Number: 98025-1-RR



Basic Information

Catalog Number:

98025-1-RR

Size:

100ug , 1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_019388

GeneID (NCBI):

12524

UNIPROT ID:

P42082-1

Full Name:

CD86 antigen

Calculated MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

230476B8

Applications

Tested Applications:

FC

Species Specificity:

mouse

Storage

Storage:

Store at 2-8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



1x10⁶ unstimulated or LPS stimulated mouse splenocytes were surface stained with 0.25 ug Anti-Mouse CD86 Rabbit Recombinant Antibody (98025-1-RR, Clone:230476B8) and PE-Conjugated AffiniPure Goat Anti-Rabbit IgG (H+L). Cells were not fixed.



1x10⁶ unstimulated (black) or LPS stimulated mouse splenocytes were surface stained with 0.25 ug Anti-Mouse CD86 Rabbit Recombinant Antibody (98025-1-RR, Clone:230476B8) and PE-Conjugated AffiniPure Goat Anti-Rabbit IgG (H+L) (red), or 0.25 ug Isotype Control (blue). Cells were not fixed.



Biolayer interferometry (BLI) kinetic assays of 98025-1-RR against Mouse CD86 were performed. The affinity constant is 25.4 nM.