

For Research Use Only

# Anti-Human CD97 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98094-2-PBS



## Basic Information

Catalog Number:

98094-2-PBS

Size:

1mg, 2 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC026690

GeneID (NCBI):

976

UNIPROT ID:

P48960

Full Name:

CD97 molecule

Calculated MW:

722 aa, 79 kDa

Purification Method:

Protein A purification

CloneNo.:

241335D8

## Applications

Tested Applications:

FC

Species Specificity:

human

## Background Information

CD97 is a member of the seven-trans-membrane subfamily of the class B G protein-coupled receptor (GPCR) group of the epidermal growth factor (EGF) and is present on the surface of lymphocytes, monocytes, macrophages, dendritic cells, granulocytes, and smooth muscle, being a dimeric glycoprotein with a 75-90 kDa intracellular domain and a 28 kDa extracellular domain (PMID: 23969601). It plays a vital role in cell adhesion, migration, and cell connection regulation (PMID: 36464675). CD97 is expressed on T-cells, but rarely on B cells. CD97 might play key roles in cellular adhesion, through connections with other proteins of the cell surface and of the extracellular matrix (PMID: 32781778).

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> human PBMCs were surface stained with 0.15 ug Anti-Human CD97 Rabbit Recombinant Antibody (98094-2-RR, Clone: 241335D8) or Isotype Control and PE-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). Cells were co-stained with APC Anti-Human CD3 (OKT3) Mouse IgG2a Recombinant Antibody (APC-65569, Clone: OKT3). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed. This data was developed using the same



1x10<sup>6</sup> human PBMCs were surface stained with 0.15 ug Anti-Human CD97 Rabbit Recombinant Antibody (98094-2-RR, Clone: 241335D8) or Isotype Control and PE-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed. This data was developed using the same antibody clone with 98094-2-PBS in a different storage buffer formulation.