

For Research Use Only

Anti-Human IL3RA/CD123 Rabbit Recombinant Antibody

Catalog Number: 98107-1-RR



Basic Information

Catalog Number:	98107-1-RR	GenBank Accession Number:	BC035407	Purification Method:	Protein A purification
Size:	100ug, 1000 µg/ml	GeneID (NCBI):	3563	CloneNo.:	241268E3
Source:	Rabbit	ENSEMBL Gene ID:	ENSG00000185291		
Isotype:	IgG	UNIPROT ID:	P26951		
		Full Name:	interleukin 3 receptor, alpha (low affinity)		
		Calculated MW:	378 aa, 43 kDa		

Applications

Tested Applications:
FC
Species Specificity:
human

Background Information

Interleukin-3 receptor subunit alpha (IL3RA), also known as CD123, belongs to the cytokine receptor family. It is a glycoprotein containing an extracellular domain, involving a predicted Ig-like domain, two FnIII domains, a transmembrane domain and an intracellular domain (PMID: 24513123). IL3RA is expressed on hematopoietic progenitor cells, plasmacytoid dendritic cells, basophils, monocytes, some B cells, and in a variety of hematologic malignancies (PMID: 10527461; 34855461). IL-3 initially binds to IL3RA and subsequently recruits the beta chain (CD131) to form the high-affinity receptor, resulting in the activation of a series of signaling pathways including the JAK/STAT, Ras-MAPK, and phosphatidylinositol 3-kinase pathways (PMID: 29296749). IL-3R signaling regulates the proliferation, survival, and differentiation of hematopoietic cells.


Storage

Storage:
Store at 2 - 8°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide, pH 7.3.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



1x10⁶ human immature monocyte-derived dendritic cells were surface stained with 0.25 ug Anti-Human IL3RA/CD123 Rabbit Recombinant Antibody (98107-1-RR, Clone:241268E3) or 0.25 ug Rabbit IgG Isotype Control Recombinant Antibody (98136-1-RR, Clone: 240953C9), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were co-stained with CoraLite® Plus 647 Anti-Human CD45RO (UCHL1) (CL647-65150, Clone: UCHL1). Cells were incubated with FC



1x10⁶ human immature monocyte-derived dendritic cells were surface stained with 0.25 ug Anti-Human IL3RA/CD123 Rabbit Recombinant Antibody (98107-1-RR, Clone:241268E3) (red) or 0.25 ug Rabbit IgG Isotype Control Recombinant Antibody (98136-1-RR, Clone: 240953C9) (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.



Biolayer interferometry (BLI) kinetic assays of 98107-1-RR against Human IL3RA/CD123 were performed. The affinity constant is 2.63 nM.