

For Research Use Only

Anti-Human TNFRSF9/CD137 Rabbit Recombinant Antibody

Catalog Number: 98112-1-RR



Basic Information

Catalog Number:

98112-1-RR

Size:

100ug , 1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC006196

GeneID (NCBI):

3604

UNIPROT ID:

Q07011

Full Name:

tumor necrosis factor receptor
superfamily, member 9

Calculated MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

241344F3

Applications

Tested Applications:

FC

Species Specificity:

human

Background Information

CD137, also known as TNFRSF9 or 4-1BB, is an inducible T cell surface receptor that belongs to the tumor necrosis factor receptor superfamily. CD137 is a transmembrane protein expressed on the surface of activated T-cells. In addition, activation-dependent expression of CD137 has also been found in B lymphocytes, monocytes, and diverse nonlymphoid cell types. CD137 provides a co-stimulatory signal that enhances the survival, and differentiation of cells, and has a crucial role in the development of CD8 cytotoxic T cells and anti-tumor immunity. (PMID: 9826581; 23696891)

Storage

Storage:

Store at 2 - 8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH 7.3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.



Selected Validation Data




1x10⁶ unstimulated or Con-A-stimulated human PBMCs were surface stained with 0.25 ug Anti-Human TNFRSF9/CD137 Rabbit Recombinant Antibody (98112-1-RR, Clone: 241344F3) and PE-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L), and 0.25 ug APC Anti-Human CD8 (UCHT4) Mouse IgG2a Recombinant Antibody (APC-65559, Clone: UCHT4). Cells were not fixed.



1x10⁶ unstimulated (black) or Con-A-stimulated human PBMCs (red) were surface stained with 0.25 ug Anti-Human TNFRSF9/CD137 Rabbit Recombinant Antibody (98112-1-RR, Clone: 241344F3) and PE-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L), or Isotype Control (blue). Cells were not fixed.



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded human tonsillitis tissue using TNFRSF9/CD137 antibody (98112-1-RR, Clone: 241344F3) at dilution of 1:180 and Multi-rAb CoraLite® Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Biolayer interferometry (BLI) kinetic assays of 98112-1-RR against Human TNFRSF9/CD137 were performed. The affinity constant is 10.3 nM.