

For Research Use Only

Anti-Mouse P-Selectin/CD62P Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98302-2-PBS



Basic Information

Catalog Number:

98302-2-PBS

Size:

1mg, 2 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_011347.2

GeneID (NCBI):

20344

UNIPROT ID:

Q01102

Full Name:

selectin, platelet

Calculated MW:

83 kDa

Purification Method:

Protein A purification

CloneNo.:

240766F5

Applications

Tested Applications:

FC

Species Specificity:

mouse

Background Information

P-selectin, also known as GMP140 or CD62P, is a transmembrane glycoprotein that mediates the interaction of activated endothelial cells or platelets with leukocytes. It is an adhesion molecule involved in the pathogenesis of inflammation, thrombosis, and oncogenesis. P-selectin is stored in the alpha-granules of platelets and Weibel-Palade bodies of endothelial cells. Upon cell activation by agonists, P-selectin is transported rapidly to the cell surface.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ Thrombin treated mouse peripheral blood platelets were surface stained with 0.25 ug Anti-Mouse P-Selectin/CD62P Rabbit RecAb (98302-2-RR, Clone: 240766F5) (red) or Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue), and PE-conjugated Goat Anti-Rabbit IgG. Untreated mouse peripheral blood platelets were surface stained with 0.25 ug Anti-Mouse P-Selectin/CD62P Rabbit RecAb (98302-2-RR, Clone: 240766F5) and PE-conjugated Goat

