

For Research Use Only

Anti-Human EPCAM/CD326 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98303-5-PBS



Basic Information

Catalog Number: 98303-5-PBS	GenBank Accession Number: NM_002354.3	Purification Method: Protein A purification
Size: 1mg, 3 mg/ml	GeneID (NCBI): 4072	CloneNo.: 250626C6
Source: Rabbit	UNIPROT ID: P16422	
Isotype: IgG	Full Name: epithelial cell adhesion molecule	
Immunogen Catalog Number: EG1370	Calculated MW: 35 kDa	

Applications

Tested Applications:
FC

Species Specificity:
human

Background Information

Epithelial cell adhesion molecule (EpCAM, CD326) is a type I transmembrane glycoprotein that functions as a homophilic, epithelial-specific intercellular cell-adhesion molecule. In addition to cell adhesion, EpCAM is also involved in cellular signaling, cell migration, proliferation, and differentiation. EpCAM is highly expressed on most carcinomas and therefore of potential use as a diagnostic and prognostic marker for a variety of carcinomas, and has become a therapeutic target. EpCAM may occur in distinct forms due to glycosylation. (PMID: 20837599; 19249674; 21576002; 22647938; 12691820)


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ HT-29 (left) or HeLa (right) cells were surface stained with 0.25 ug Anti-Human EPCAM/CD326 Rabbit RecAb (98303-5-RR, Clone: 250626C6) (red) or Rabbit IgG Isotype Control Recombinant Antibody (98136-1-RR, Clone: 240953C9) (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed. This data was developed using the same antibody clone with 98303-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 98303-5-RR against Human EPCAM/CD326 were performed. The affinity constant is 14.4 nM.