

For Research Use Only

Anti-Mouse CD164 Rabbit Recombinant Antibody

Catalog Number: 98419-3-RR



Basic Information

Catalog Number:

98419-3-RR

Size:

100ug, 1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2645

GenBank Accession Number:

NM_016898.2

GeneID (NCBI):

53599

UNIPROT ID:

Q9ROL9

Full Name:

CD164 antigen

Calculated MW:

21 kDa

Purification Method:

Protein A purification

CloneNo.:

242986C2

Recommended Dilutions:

FC: 0.25 ug per 10⁶ cells in a 100 µl suspension

Applications

Tested Applications:

FC

Species Specificity:

mouse

Positive Controls:

FC : Transfected HEK-293T cells,

Storage

Storage:

Store at 2 - 8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ HEK-293T cells or CD164 transfected HEK-293T cells were surface stained with 0.25 ug Anti-Mouse CD164 Rabbit RecAb (98419-3-RR, Clone:242986C2) and Multi-rAb CoraLite® Plus 647-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR005). Cells were not fixed.

1x10⁶ HEK-293T cells or CD164-transfected HEK-293T cells were surface stained with 0.25 ug Anti-Mouse CD164 Rabbit RecAb (98419-3-RR, Clone: 242986C2) (red) or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue), and Multi-rAb CoraLite® Plus 647-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR005). Cells were not fixed.