

For Research Use Only

Anti-Rabbit IgM Mouse Recombinant Antibody

Catalog Number: 98426-2-MR



Basic Information

Catalog Number:

98426-2-MR

Size:

100ug , 1000 µg/ml

Source:

Mouse

Isotype:

IgG

Immunogen Catalog Number:

EG4359

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

242572B6

Recommended Dilutions:

FC: 0.25 ug per 10⁶ cells in a 100 µl suspension

Applications

Tested Applications:

FC

Species Specificity:

rabbit

Positive Controls:

FC : rabbit splenocytes,

Storage

Storage:

Store at 2 - 8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ rabbit splenocytes were surface stained with 0.25 ug Anti-Rabbit IgM Mouse RecAb (98426-2-MR, Clone: 242572B6) or Mouse IgG2a Isotype Control (C1.18.4) (65208-1-Ig, Clone: C1.18.4), and PE-conjugated Goat Anti-Mouse IgG. Cells were then stained with CoraLite® Plus 488 Anti-Rabbit CD3. Cells were not fixed.



1x10⁶ rabbit splenocytes were surface stained with 0.25 ug Anti-Rabbit IgM Mouse RecAb (98426-2-MR, Clone: 242572B6) (red) or Mouse IgG2a Isotype Control (C1.18.4) (65208-1-Ig, Clone: C1.18.4) (blue), and PE-conjugated Goat Anti-Mouse IgG. Cells were then stained with CoraLite® Plus 488 Anti-Rabbit CD3. Cells were not fixed.