

For Research Use Only

Anti-Mouse CD4 Rabbit Recombinant Antibody

Catalog Number: 98522-3-RR



Basic Information

Catalog Number: 98522-3-RR	GenBank Accession Number: NM_013488	Purification Method: Protein A purification
Size: 100ug , 1000 µg/ml	GeneID (NCBI): 12504	CloneNo.: 250263B7
Source: Rabbit	UNIPROT ID: P06332	Recommended Dilutions: FC: 0.125 ug per 10 ⁶ cells in a 100 µl suspension
Isotype: IgG	Full Name: CD4 antigen	
Immunogen Catalog Number: EG3847	Calculated MW: 51 kDa	

Applications

Tested Applications: FC	Positive Controls: FC : mouse splenocytes,
Species Specificity: mouse	

Background Information

CD4 is a 55-kDa transmembrane glycoprotein expressed on T helper cells, majority of thymocytes, monocytes, macrophages, and dendritic cells (PMID: 9304802; 12213222). CD4 is an accessory protein for MHC class-II antigen/T-cell receptor interaction. It plays an important role in T helper cell development and activation (PMID: 9539765; 3112582). CD4 serves as a receptor for the human immunodeficiency virus (HIV) (PMID: 9304802).


Storage

Storage:
Store at 2 - 8°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



1x10⁶ mouse splenocytes were surface stained with 0.125 ug Anti-Mouse CD4 Rabbit RecAb (98522-3-RR, Clone: 250263B7) or Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were then stained with CoraLite® Plus 647 Anti-Mouse CD3 (17A2) (CL647-65077, Clone: 17A2). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.



1x10⁶ mouse splenocytes were surface stained with 0.125 ug Anti-Mouse CD4 Rabbit RecAb (98522-3-RR, Clone: 250263B7) (red) or Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.



Biolayer interferometry (BLI) kinetic assays of 98522-3-RR against Mouse CD4 were performed. The affinity constant is below 1 pM.