

For Research Use Only

Anti-Mouse IL-25/IL-17E Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98536-3-PBS



Basic Information

Catalog Number: 98536-3-PBS	GenBank Accession Number: NM_080729.3	Purification Method: Protein A purification
Size: 1mg, 3 mg/ml	GeneID (NCBI): 140806	CloneNo.: 250033D3
Source: Rabbit	UNIPROT ID: Q8VHH8	
Isotype: IgG	Full Name: interleukin 25	
Immunogen Catalog Number: EG3316	Calculated MW: 19KD	

Applications

Tested Applications:
FC (Intra)

Species Specificity:
mouse

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ HEK-293T cells or IL-25/IL-17E transfected HEK-293T cells were intracellularly stained with 0.25 ug Anti-Mouse IL-25/IL-17E Rabbit RecAb (98536-3-RR, Clone: 250033D3) and Multi-rAb CoraLite[®] Plus 647-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR005). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer. This data was developed using the same antibody clone with 98536-3-PBS in a different storage



1x10⁶ HEK-293T cells or IL-25/IL-17E transfected HEK-293T cells were intracellularly stained with 0.25 ug Anti-Mouse IL-25/IL-17E Rabbit RecAb (98536-3-RR, Clone: 250033D3) (red) or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue), and Multi-rAb CoraLite[®] Plus 647-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR005). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry

