

For Research Use Only

Anti-Mouse IL-25/IL-17E Rabbit Recombinant Antibody

Catalog Number: 98536-3-RR



Basic Information

Catalog Number:

98536-3-RR

Size:

100ug, 1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3316

GenBank Accession Number:

NM_080729.3

GeneID (NCBI):

140806

UNIPROT ID:

Q8VHH8

Full Name:

interleukin 25

Calculated MW:

19KD

Purification Method:

Protein A purification

CloneNo.:

250033D3

Recommended Dilutions:

FC (Intra): 0.25 ug per 10⁶ cells in a 100 µl suspension

Applications

Tested Applications:

FC (Intra)

Species Specificity:

mouse

Positive Controls:

FC (Intra) : Transfected HEK-293T cells,

Storage

Storage:

Store at 2 - 8°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ HEK-293T cells or IL-25/IL-17E transfected HEK-293T cells were intracellularly stained with 0.25 ug Anti-Mouse IL-25/IL-17E Rabbit RecAb (98536-3-RR, Clone: 250033D3) and Multi-rAb CoraLite[®] Plus 647-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR005). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer.

1x10⁶ HEK-293T cells or IL-25/IL-17E transfected HEK-293T cells were intracellularly stained with 0.25 ug Anti-Mouse IL-25/IL-17E Rabbit RecAb (98536-3-RR, Clone: 250033D3) (red) or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue), and Multi-rAb CoraLite[®] Plus 647-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR005). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer.