

For Research Use Only

Anti-Human CD21 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98664-5-PBS



Basic Information

Catalog Number:	98664-5-PBS	GenBank Accession Number:	NM_001877.5	Purification Method:	Protein A purification
Size:	1mg, 3 mg/ml	GeneID (NCBI):	1380	CloneNo.:	252181G2
Source:	Rabbit	ENSEMBL Gene ID:	ENSG00000117322		
Isotype:	IgG	UNIPROT ID:	P20023-1		
Immunogen Catalog Number:	EG2980	Full Name:	complement component (3d/Epstein Barr virus) receptor 2		
		Calculated MW:	113 kDa		

Applications

Tested Applications:
FC
Species Specificity:
human

Background Information

CD21, also known as complement receptor type 2 (CR2), complement C3d receptor, and Epstein-Barr virus receptor, is a transmembrane protein that contains a small cytoplasmic domain, a transmembrane region and an extracellular domain consisting of 15 tandem short consensus repeat sequences (PMID: 6230668; 2551147). It is expressed on B cells, follicular dendritic cells, thymocytes and a subset of peripheral T cells (PMID: 26119182). CD21 binds complement fragments C3d, C3dg and iC3b and acts as a receptor for the Epstein-Barr virus (PMID: 7753047). On B cells, together with CD19 and CD81, CD21 forms a complex that functions as a co-receptor to the BCR (PMID: 7542009). It is involved in B cells activation and functions.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



1x10⁶ human PBMCs were surface stained with 0.25 ug Anti-Human CD21 Rabbit RecAb (98664-5-RR, Clone: 252181G2) or 0.25 ug Rabbit IgG Isotype Control Recombinant Antibody (98136-1-RR, Clone: 240953C9), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were then stained with CoraLite® Plus 647 Anti-Human CD19 (4G7) Mouse IgG2a Recombinant Antibody (CL647-65562, Clone: 4G7). Cells were not fixed. Lymphocytes were gated. This data was



1x10⁶ human PBMCs were surface stained with 0.25 ug Anti-Human CD21 Rabbit RecAb (98664-5-RR, Clone: 252181G2) (red) or 0.25 ug Rabbit IgG Isotype Control Recombinant Antibody (98136-1-RR, Clone: 240953C9) (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed. Lymphocytes were gated. This data was developed using the same antibody clone with 98664-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 98664-5-RR against Human CD21 were performed. The affinity constant is 35.0 pM.