

For Research Use Only

# Anti-Human Claudin 19 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98681-1-PBS



## Basic Information

Catalog Number:

98681-1-PBS

Size:

1mg, 3 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC030524

GeneID (NCBI):

149461

UNIPROT ID:

Q8N6F1

Full Name:

claudin 19

Calculated MW:

23 kDa

Purification Method:

Protein A purification

CloneNo.:

251126D7

## Applications

Tested Applications:

FC

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> U2OS cells were surface stained with 0.25 ug Anti-Human Claudin 19 Rabbit RecAb (98681-1-RR, Clone: 251126D7) (red) or 0.25 ug Isotype Control (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed. This data was developed using the same antibody clone with 98681-1-PBS in a different storage buffer formulation.