

For Research Use Only

Mouse IgG1, λ isotype control Recombinant Antibody, PBS Only

Catalog Number: 98741-1-PBS



Basic Information

Catalog Number:

98741-1-PBS

Size:

1mg, 3 mg/ml

Source:

Mouse

Isotype:

IgG1, lambda

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

253505A1

Applications

Tested Applications:

FC

Species Specificity:

n/a

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with 0.25 ug Mouse IgG1, λ isotype control RecAb (98741-1-MR, Clone:253505A1) (red in histogram) and PE-conjugated Goat Anti-Mouse IgG, or 0.25 ug Anti-Human CD4 (RPA-T4) (65143-1-Ig, Clone: RPA-T4) (green in histogram), or unstained cells (blue in histogram). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed. This data was developed using the same antibody clone with 98741-1-PBS in a



1x10⁶ rat splenocytes were surface stained with 0.25 ug Mouse IgG1, λ isotype control RecAb (98741-1-MR, Clone:253505A1) (red in histogram) and PE-conjugated Goat Anti-Mouse IgG, or 0.25 ug Anti-Rat CD45 (OX1) (65228-1-Ig, Clone: OX1) (green in histogram), or unstained cells (blue in histogram). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed. This data was developed using the same antibody clone with 98741-1-PBS in a

