

For Research Use Only

Anti-Rat CCL4/MIP-1 beta Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98770-4-PBS



Basic Information

Catalog Number:

98770-4-PBS

Size:

1mg, 3 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3200

GenBank Accession Number:

NM_053858.1

GeneID (NCBI):

116637

UNIPROT ID:

P50230

Full Name:

chemokine (C-C motif) ligand 4

Calculated MW:

10 kDa

Purification Method:

Protein A purification

CloneNo.:

251973F1

Applications

Tested Applications:

FC (Intra)

Species Specificity:

rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3



For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ untreated or LPS-treated NR8383 cells were intracellularly stained with 0.25 ug Anti-Rat CCL4/MIP-1 beta Rabbit RecAb (98770-4-RR, Clone: 251973F1) and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were fixed and permeabilized with Intracellular Flow Cytometry Fixation & Permeabilization Buffer Kit (PF00019). This data was developed using the same antibody clone with 98770-4-PBS in a different storage buffer formulation.

1x10⁶ LPS treated NR8383 cells were intracellularly stained with 0.25 ug Anti-Rat CCL4/MIP-1 beta Rabbit RecAb (98770-4-RR, Clone: 251973F1)(red) or Rabbit IgG Isotype Control Recombinant Antibody (98136-1-RR, Clone: 240953C9) (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). 1x10⁶ untreated NR8383 were surface stained with 0.25 ug Anti-Rat CCL4/MIP-1 beta Rabbit RecAb (98770-4-RR, Clone: 251973F1) and PE-Conjugated Goat Anti-Rabbit IgG(H+L)

