

For Research Use Only

Anti-Human uPA/Urokinase Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98824-2-PBS



Basic Information

Catalog Number:

98824-2-PBS

Size:

1mg, 3 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG7083

GenBank Accession Number:

NM_002658.6

GeneID (NCBI):

5328

UNIPROT ID:

P00749-1

Full Name:

plasminogen activator, urokinase

Calculated MW:

49kDa

Purification Method:

Protein A purification

CloneNo.:

252624A9

Applications

Tested Applications:

FC

Species Specificity:

human

Background Information

uPA (Urokinase-type plasminogen activator) is also named as Urokinase or PLAU, and belongs to the peptidase S1 family. uPA initiates a proteolytic cascade, which degrades the extracellular matrix during tissue growth and remodelling. The secreted proform of uPA is activated after binding to a specific high-affinity cell surface receptor. Secreted UPA is exposed to rapid inactivation in the pericellular space by plasminogen activator inhibitor-1 (PAI-1), which is stored in the extracellular matrix bound to vitronectin (PMID: 9665342). It has 2 isoforms produced by alternative splicing. Defects in uPA are the cause of Quebec platelet disorder (QPD).

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ HT-1080 (left) or Ramos cells (right) were surface stained with 0.25 ug Anti-Human uPA/Urokinase Rabbit RecAb (98824-2-RR, Clone: 252624A9) (red) or 0.25 ug Rabbit IgG Isotype Control Recombinant Antibody (98136-1-RR, Clone: 240953C9) (blue), and PE-conjugated Goat Anti-Rabbit IgG. Cells were not fixed. This data was developed using the same antibody clone with 98824-2-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 98824-2-RR against Human uPA/Urokinase were performed. The affinity constant is below 1 pM.