

For Research Use Only

Atlantic Blue™ Anti-Mouse LAG-3/CD223 (C9B7W)



Catalog Number: AB-65098

Basic Information

Catalog Number:

AB-65098

Size:

100ug , 500 µg/ml

Source:

Rat

Isotype:

IgG1, kappa

GenBank Accession Number:

BC120591

GeneID (NCBI):

16768

Full Name:

Lymphocyte-activation gene 3

Purification Method:

Affinity purification

CloneNo.:

C9B7W

Excitation/Emission maxima wavelengths:

404 nm / 458 nm

Applications

Tested Applications:

FC

Species Specificity:

mouse

Background Information

Lymphocyte activation gene 3 protein (LAG-3), also known as CD223, is a type I transmembrane protein belonging to the immunoglobulin superfamily (PMID: 1692078). It is a CD4-related molecule that binds MHC class II with high affinity (PMID: 1692078; 9159144). LAG-3 is expressed on activated T and NK cells. Besides, natural CD4(+)CD25(+) Tregs express LAG-3 upon activation (PMID: 15485628). LAG-3 negatively regulates T-cell function and homeostasis, and contributes to the suppressor activity of Tregs (PMID: 14644131; 15485628).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ plate-bound Anti-CD3 and Anti-CD28 stimulated mouse splenocytes were surface co-stained with CoraLite®488 Plus Anti-Mouse CD3 (17A2) (CL488-65077, Clone: 17A2) and 0.5 ug Atlantic Blue™ Anti-Mouse CD223 (AB-65098, Clone: C9B7W), or 0.5 ug Isotype Control. Cells were not fixed.