

For Research Use Only

Atlantic Blue™ Anti-Human CD1a (HI149)

Catalog Number: AB-65172



Basic Information

Catalog Number:

AB-65172

Size:

100 tests, 5 ul/test

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

BC031645

GeneID (NCBI):

909

UNIPROT ID:

P06126

Full Name:

CD1a molecule

Calculated MW:

327 aa, 37 kDa

Purification Method:

Unknown

CloneNo.:

HI149

Excitation/Emission maxima
wavelengths:

404 nm / 458 nm

Applications

Tested Applications:

FC

Species Specificity:

human

Background Information

CD1a is an antigen-presenting protein that binds self and non-self lipid and glycolipid antigens and presents them to T-cell receptors on natural killer T-cells. During protein synthesis and maturation, CD1 family members bind endogenous lipids that are replaced by lipid or glycolipid antigens when the proteins are internalized and pass through endosomes, before trafficking back to the cell surface.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ MOLT-4 cells were surface stained with 5 ul Atlantic Blue™ Anti-Human CD1a (AB-65172, Clone: HI149) or Atlantic Blue™ Mouse IgG1 Isotype Control (AB-65124, Clone: MOPC-21). Cells were not fixed.



1x10⁶ MOLT-4 cells were surface stained with 5 ul Atlantic Blue™ Anti-Human CD1a (AB-65172, Clone: HI149) or Atlantic Blue™ Mouse IgG1 Isotype Control (AB-65124, Clone: MOPC-21). Cells were not fixed.