

For Research Use Only

Atlantic Blue™ Anti-Human CD146 (P1H12)



Catalog Number: **AB-65181**

Basic Information

Catalog Number: AB-65181	GenBank Accession Number: BC056418	Purification Method: Affinity purification
Size: 100tests , 5 ul/test	GeneID (NCBI): 4162	CloneNo.: P1H12
Source: Mouse	Full Name: melanoma cell adhesion molecule	Excitation/Emission maxima wavelengths: 404 nm / 458 nm
Isotype: IgG1, kappa	Calculated MW: 646 aa, 72 kDa	

Applications

Tested Applications:
FC

Species Specificity:
Rabbit, Canine, Mouse, Human

Background Information

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ A375 cells were surface stained with 5 ul Atlantic Blue™ Anti-Human CD146 (AB-65181, Clone:P1H12), or Mouse IgG1 Isotype Control. Cells were not fixed.



1X10⁶ A375 cells were surface stained with 5 ul Atlantic Blue™ Anti-Human CD146 (AB-65181, Clone:P1H12), or Mouse IgG1 Isotype Control. Cells were not fixed.