

For Research Use Only

Atlantic Blue™ Anti-Rat CD43 (W3/13)



Catalog Number: **AB-65237**

Basic Information

Catalog Number: AB-65237	GenBank Accession Number: GeneID (NCBI): 24796	Purification Method: Affinity purification
Size: 100ug, 0.5 mg/ml	Full Name: sialophorin	CloneNo.: W3/13
Source: Mouse	Calculated MW: 38 kDa	Excitation/Emission maxima wavelengths: 404 nm / 458 nm
Isotype: IgG1, kappa		

Applications

Tested Applications:
FC
Species Specificity:
Rat

Background Information


Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ Wistar rat splenocytes were surface stained with 0.5 ug Atlantic Blue™ Anti-Rat CD43 (AB-65237, Clone: W3/13) or 0.5 ug Control Antibody. Cells were not fixed.