

For Research Use Only

# Atlantic Blue™ Anti-Mouse CD94 (18D3)



Catalog Number: **AB-65240**

## Basic Information

Catalog Number:

AB-65240

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC120854

GeneID (NCBI):

16643

UNIPROT ID:

O54707

Full Name:

killer cell lectin-like receptor,  
subfamily D, member 1

Purification Method:

Affinity purification

CloneNo.:

18D3

Excitation/Emission maxima  
wavelengths:

404 nm / 458 nm

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> mouse splenocytes were surface co-stained with FITC Plus Anti-Mouse CD49b (HMA2) (FITC-65245, Clone: HMA2) and 0.5 ug Atlantic Blue™ Anti-Mouse CD94 (AB-65240, Clone: 18D3), or 0.5 ug Isotype Control. Cells were not fixed.