

For Research Use Only

# Atlantic Blue™ Mouse IgM Isotype Control (MOPC-104E)



Catalog Number: **AB-65243**

## Basic Information

Catalog Number: AB-65243	GenBank Accession Number: GeneID (NCBI):	Purification Method: N/A
Size: 100ug, 200 µg/ml	Full Name:	CloneNo.: MOPC-104E
Source: Mouse		Excitation/Emission maxima wavelengths: 404 nm / 458 nm
Isotype: IgM, lambda		

## Applications

Tested Applications:  
FC  
Species Specificity:

## Background Information


### Storage

Storage:  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> human whole blood cells were surface stained with 0.2 ug Atlantic Blue™ Anti-Mouse IgM Isotype Control (red), 5 ul Atlantic Blue™ Anti-Human STRO-1 (green), or unstained (blue). Cells were not fixed. CD45dim, CD34+ lymphocytes were gated.