

For Research Use Only

Atlantic Blue™ Anti-Mouse CD49b (HMA2)



Catalog Number: AB-65245

1 Publications

Basic Information

Catalog Number:

AB-65245

Size:

100ug , 500 µg/ml

Source:

Armenian Hamster

Isotype:

IgG

GenBank Accession Number:

BC065139

GeneID (NCBI):

16398

UNIPROT ID:

Q62469

Full Name:

integrin alpha 2

Purification Method:

Affinity purification

CloneNo.:

HMA2

Excitation/Emission maxima wavelengths:

404 nm / 458 nm

Applications

Tested Applications:

FC

Cited Applications:

FC

Species Specificity:

Mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Dunfang Wang	39034016	J Ethnopharmacol	FC

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface stained with 0.5 ug Atlantic Blue™ Anti-Mouse CD49b (AB-65245, Clone:HMA2) or 0.5 ug Control Antibody. Cells were not fixed.