

For Research Use Only

APC Anti-Mouse CD3 ϵ (500-A2)

Catalog Number: APC-65061



Basic Information

Catalog Number: APC-65061	GenBank Accession Number: BC098236	Purification Method: Affinity purification
Size: 100ug, 0.1 mg/ml	GeneID (NCBI): 12501	CloneNo.: 500-A2
Source: Hamster	Full Name: CD3 antigen, epsilon polypeptide	Excitation/Emission maxima wavelengths: 650 nm / 660 nm
Isotype: IgG		

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

CD3 is a multimeric protein associated with the T-cell receptor (TCR) to form a complex involved in antigen recognition and signal transduction (PMID: 15885124). CD3 is composed of CD3 γ , δ , ϵ , and ζ chains (PMID: 1826255). It is expressed by thymocytes in a developmentally regulated manner, T cells, and some NK cells (PMID: 3289580). The TCR recognizes antigens bound to major histocompatibility complex (MHC) molecules. TCR-mediated peptide-MHC recognition is transmitted to the CD3 complex, leading to the intracellular signal transduction (PMID: 11985657).

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and a stabilizer.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 C57BL/6 mouse splenocytes were surface stained with 0.20 μ g APC-Anti-Mouse CD3 ϵ (APC-65061, clone 500-A2) (red) or 0.20 μ g APC-Hamster IgG isotype control (blue). Cells were not fixed.