

For Research Use Only

# APC Anti-Mouse ICAM-1/CD54 (YN1/1.7.4)

Catalog Number: APC-65066

1 Publications



## Basic Information

Catalog Number:

APC-65066

Size:

100ug, 0.1 mg/ml

Source:

Rat

Isotype:

IgG2b, kappa

GenBank Accession Number:

BC008626

GeneID (NCBI):

15894

UNIPROT ID:

P13597

Full Name:

intercellular adhesion molecule 1

Purification Method:

Affinity purification

CloneNo.:

YN1/1.7.4

Excitation/Emission maxima  
wavelengths:

650 nm / 660 nm

## Applications

Tested Applications:

FC

Cited Applications:

FC

Species Specificity:

mouse

Cited Species:

mouse

## Background Information

ICAM-1 (CD54) is a 90-kDa transmembrane glycoprotein of the immunoglobulin superfamily and is critical for the firm attachment and transmigration of leukocytes out of blood vessels and into tissues (PMID: 19307690). ICAM-1 is expressed by several cell types, typically on endothelial cells and cells of the immune system, and its expression can be up-regulated by various stimuli, including TNF- $\alpha$ , INF- $\gamma$ , IL-1 and thrombin (PMID: 3086451; 9694714; 15979056). It is a ligand for LFA-1 and Mac-1, serves as a receptor for rhinovirus, and is one of several receptors used by Plasmodium falciparum (PMID: 2566624; 2538244; 2475784).

## Notable Publications

Author	Pubmed ID	Journal	Application
Zhenzhen Zhang	34270853	Small	FC

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and a stabilizer.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> C57BL/6 mouse splenocytes were surface stained with 0.20 ug APC-Anti-Mouse CD54 (APC-65066, clone YN1/1.7.4) (red) or 0.20 ug APC-rat IgG2b isotype control (blue). Cells were not fixed.