

For Research Use Only

# APC Anti-Human CD1a (HI149)

Catalog Number: APC-65172



## Basic Information

Catalog Number:

APC-65172

Size:

100 tests, 10 µl/test

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

BC031645

GeneID (NCBI):

909

UNIPROT ID:

P06126

Full Name:

CD1a molecule

Calculated MW:

327 aa, 37 kDa

Purification Method:

The purified antibody is conjugated with allophycocyanin (APC) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

CloneNo.:

HI149

Excitation/Emission maxima wavelengths:

650 nm / 660 nm

## Applications

Tested Applications:

FC

Species Specificity:

Human

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> MOLT-4 cells were surface stained with 10 ul APC Anti-Human CD1a (APC-65172, Clone:HI149), or Mouse IgG1 Isotype Control. Cells were not fixed.



1X10<sup>6</sup> MOLT-4 cells were surface stained with 10 ul APC Anti-Human CD1a (APC-65172, Clone:HI149), or Mouse IgG1 Isotype Control. Cells were not fixed.