

For Research Use Only

APC Anti-Human IgD (IA6-2)

Catalog Number: **APC-65175**



Basic Information

Catalog Number:

APC-65175

Size:

100tests , 10 µl/test

Source:

Mouse

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC021276

GeneID (NCBI):

3495

Full Name:

immunoglobulin heavy constant
delta

Calculated MW:

573 aa, 63 kDa

Purification Method:

The purified antibody is conjugated with allophycocyanin (APC) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

CloneNo.:

IA6-2

Excitation/Emission maxima wavelengths:

650 nm / 660 nm

Applications

Tested Applications:

FC

Species Specificity:

Human

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ human peripheral blood lymphocytes were surface stained with PE Anti-Human CD19 (PE-65110, Clone:H1B19) and 10.00 uL/test APC Anti-Human IgD (APC-65175, Clone:IA6-2). Sample were not fixed.



1X10⁶ human peripheral blood lymphocytes were surface stained with PE Anti-Human CD19 (PE-65110, Clone:H1B19) and APC mouse IgG2a isotype control antibody. Sample were not fixed.