

For Research Use Only

APC Anti-Human TIGIT Rabbit Recombinant Antibody

Catalog Number: APC-98023-2



Basic Information

Catalog Number:

APC-98023-2

Size:

100 tests, 5 ul/test

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0749

GenBank Accession Number:

NM_173799.4

GeneID (NCBI):

201633

UNIPROT ID:

Q495A1-1

Full Name:

T cell immunoreceptor with Ig and ITIM domains

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

240423C7

Recommended Dilutions:

FC: 5 ul per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

650 nm / 660 nm

Applications

Tested Applications:

FC

Species Specificity:

human

Positive Controls:

FC : human PBMCs,

Background Information

TIGIT (T-cell immunoreceptor with Ig and ITIM domains), also known as VSIG9 or VSTM3, is an immune receptor expressed on T cells, including Treg and memory subsets, as well as on NK cells (PMID: 19011627). It contains an immunoglobulin variable domain, a transmembrane domain and an immunoreceptor tyrosine-based inhibitory motif (ITIM). TIGIT binds to poliovirus receptor (PVR, also called CD155) with high affinity, and also to PVRL2 (CD112) with lower affinity. The interaction of TIGIT with PVR on dendritic cells increases the secretion of IL-10 and decreases the secretion of proinflammatory cytokine and suppresses T-cell activation by promoting the generation of mature immunoregulatory dendritic cells. TIGIT can inhibit NK cytotoxicity directly through its ITIM (PMID: 19815499).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with 5 ul FITC Plus Anti-Human CD45RO (UCHL1) (FITC-65150, Clone: UCHL1), and 5 ul APC Anti-Human TIGIT Rabbit RecAb (APC-98023-2, Clone: 240423C7) or APC Rabbit IgG Isotype Control RecAb (APC-98136, Clone: 240953C9). Cells were not fixed.



1x10⁶ human PBMCs were surface stained with 5 ul APC Anti-Human TIGIT Rabbit RecAb (APC-98023-2, Clone: 240423C7) or APC Rabbit IgG Isotype Control RecAb (APC-98136, Clone: 240953C9). Cells were not fixed.