

For Research Use Only

FcZero-rAb™ APC Anti-Human STRO-1 (STRO-1) Rabbit IgG Recombinant Antibody

Catalog Number: APC-FcA65585



Basic Information

Catalog Number:

APC-FcA65585

Size:

100 tests, 5 µl/test

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

STRO-1

Recommended Dilutions:

FC: 5 µl per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

650 nm / 660 nm

Applications

Tested Applications:

FC

Species Specificity:

human

Positive Controls:

FC : Human bone marrow mesenchymal stromal cells,

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human bone marrow mesenchymal stromal cells were surface stained with FcZero-rAb™ PE Anti-Human CD90, and 5 ul APC Anti-Human STRO-1 (STRO-1) Rabbit IgG RecAb (APC-FcA65585, Clone: STRO-1) or APC Rabbit IgG Isotype Control RecAb (APC-FcA98136, Clone: 240953C9). Cells were not fixed.



1x10⁶ human bone marrow mesenchymal stromal cells were surface stained with 5 ul APC Anti-Human STRO-1 (STRO-1) Rabbit IgG RecAb (APC-FcA65585, Clone: STRO-1) (red) or APC Rabbit IgG Isotype Control RecAb (APC-FcA98136, Clone: 240953C9) (blue). Cells were not fixed.