

For Research Use Only

FcZero-rAb™ APC Anti-Human TNF beta Rabbit Recombinant Antibody

Catalog Number: APC-FcA98056



Basic Information

Catalog Number:

APC-FcA98056

Size:

100 tests, 5 µl/test

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0836

GenBank Accession Number:

BC034729

GeneID (NCBI):

4049

Full Name:

Lymphotoxin alpha (TNF superfamily, FC (Intra): 5 µl per 10⁶ cells in a 100 µl suspension

Calculated MW:

205 aa, 22 kDa

Purification Method:

Protein A purification

CloneNo.:

240393C3

Recommended Dilutions:

Excitation/Emission maxima wavelengths:

650 nm / 660 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human

Positive Controls:

FC (Intra): Anti-CD3/CD28 treated human PBMCs,

Background Information

TNF beta also known as lymphotoxin alpha (LTA), is a member of the tumor necrosis factor family, and is a cytokine produced by lymphocytes. TNF Beta is highly inducible, secreted, and forms heterotrimers with lymphotoxin-beta which anchor lymphotoxin-alpha to the cell surface. TNF Beta also mediates a large variety of inflammatory, immunostimulatory, and antiviral responses, is involved in the formation of secondary lymphoid organs during development and plays a role in apoptosis.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ untreated or anti-CD3/CD28 treated human PBMCs were intracellularly stained with 5 ul APC Anti-Human TNF beta Rabbit RecAb (APC-FcA98056, Clone: 240393C3) and FcZero-rAb™ CoraLite® Plus 488 Anti-Human CD3 (UCHT1) Rabbit IgG Recombinant Antibody (CL488-FcA65570, Clone: UCHT1). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).