

For Research Use Only

# FcZero-rAb® APC Anti-Mouse Tissue Factor/CD142 Rabbit Recombinant Antibody

Catalog Number: APC-FcA98114



## Basic Information

Catalog Number:

APC-FcA98114

Size:

100ug, 200 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG1102

GenBank Accession Number:

GeneID (NCBI):

14066

UNIPROT ID:

P20352

Full Name:

coagulation factor III

Calculated MW:

33 kDa

Purification Method:

Protein A purification

CloneNo.:

240875G1

Recommended Dilutions:

FC: 0.1 ug per 10<sup>6</sup> cells in a 100 µl suspension

Excitation/Emission maximum wavelengths:

650 nm / 660 nm

## Applications

Tested Applications:

FC

Species Specificity:

mouse

Positive Controls:

FC : LPS treated RAW 264.7 cells,

## Background Information

Normal blood coagulation is a complex process, involving a cascade of activation of different plasma proteins, ultimately resulting in the formation of a clot, called fibrin. Tissue Factor (TF), also named CD142, coagulation factor III or F3, is the primary initiator of the blood coagulation cascade and plays an essential role in hemostasis (PMID: 26877187; PMID: 33796108).

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> untreated or LPS-treated RAW 264.7 cells were surface stained with 0.1 ug APC Anti-Mouse Tissue Factor/CD142 Rabbit RecAb (APC-FcA98114, Clone: 240875G1). Cells were not fixed.

1x10<sup>6</sup> LPS-treated RAW 264.7 cells were surface stained with 0.1 ug APC Anti-Mouse Tissue Factor/CD142 Rabbit RecAb (APC-FcA98114, Clone: 240875G1) (red) or 0.1 ug APC Rabbit IgG Isotype Control RecAb (APC-FcA98136, Clone: 240953C9) (blue). 1x10<sup>6</sup> untreated RAW 264.7 cells were surface stained with 0.1 ug APC Anti-Mouse Tissue Factor/CD142 Rabbit RecAb (APC-FcA98114, Clone: 240875G1) (black, dashed). Cells were not fixed.

