

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

CD320 Fusion Protein



Basic Information

Catalog Number:

Ag0344

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

SDEEECRIEPCTQKGQC P P P P L P C P C T G V S D C S G G T
D K K L R N C S R L A C L A G E L R C T L S D D C I P L T W R C D G H P D
C P D S S D E L G C G T N E I L P E G D A T T M G P P V T L E S V T S L R N
A T T M G P P V T L E S V P S V G N A T S S S A G D Q S G S P T A Y G V
I A A A A V L S A S L V T A T L L L L S W L R A Q E R L R P L G L L V A M K
E S L L L S E Q K T S

(84-280 aa encoded by BC000668)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.