

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human TMPRSS2 protein



Basic Information

Catalog Number:
Ag0435

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
MALNSGPPAIGPYENHGYQPENPYPAQPTVVPTV
YEVHPAQYYPSVPVQYAPRVLTQASNPVVCTQPKSP
SGTVCTSKTKKALCITLTLGTFLVGAALAAAGLLWKFMG
SKCSNSGIECDSSGTCINPSNWCDGVSHCPGGEDEN
RCVRLYGPNFILQVYSSQRKSWHPVCQDDWENYNG
RAACRDMGYKNNFYSSQGIVDDSGSTSFMKLNTSAG
NVDIYKKLYHSDACSSKAVVSLRCIACGVNLSRQS
RIVGGESALPGAWPWQVSLHVQNVHVC GGSII TPE
WIVTAAHCVEKPLNPNWHWTAFA GILRQSF MFYGA
GYQVEKVISHPNYDSKT KNDIALMKLQKPLTFNDLV
KPVCLPNPGMMLQPEQLCWISGWGATEEKGKTSEVL
NAAKVLIIETQRCNSRVYVDNLITPAMICAGFLQGNV
DSCQGDSSGGPLVTSKNNIWWLIGDTSWGS GCAKAY
RPGVYGNVMVFTDWIYRQMRADG
(1-492aa encoded by BC051839)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W: ptglab.com
(toll free in USA), or 1(312) 455-8498
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.