

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# TMPRSS2 Fusion Protein



## Basic Information

**Catalog Number:**  
Ag0435

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Validated Application:**  
Blocking peptide

**Peptide Sequence:**  
MALNSGPPAIGPYENHGYQPENYPQAQPTVVPTV  
YEVHPAQYYSPVPQYAPRVLTAQSNPVVCTQPKSP  
SGTVCTSKTKKALCITLTLGTFLVGAALAAAGLLWKFMG  
SKCSNSGIECDSSGTCINPSNWC DGVSHCPGGEDE  
RCVRLYGPNFILQVYSSQRKSWHPVCQDDWENY  
RAACRDMGYKNNFYSSQGIVDDSGSTSFMKLN  
TSAGNVDIYKKLYHSDACSSKAVVSLRCIACG  
VNLN SSRQSRIVGGESALPGAWPWQVSLHVQ  
NVHVCGGSIITPEWIVTAAHCEKPLNPN  
PWHWTAFAGILRQSFMFYGAGYQVEKVI  
SHPNYDSKTKNNDIALMKLQKPLTFNDLV  
KPVCLPNPGMMLQPEQLCWISGWGATEEKG  
KTSEVLNAAKVLIIETQRCSNRVYVDNLIT  
PAMICAGFLQGNVDSCQGDSSGGPLVTSK  
NNIWWLIGDTSWGS GCAKAYRPGVYGNV  
MVFTDWIYRQMRADG  
(1-492aa encoded by BC051839)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W: ptglab.com  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.