

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human NELL2 protein



## Basic Information

**Catalog Number:**  
Ag0703

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*e coli*-derived, PGEX-tag, with N-terminal N/A.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Peptide Sequence:**  
MESRVLLRTFCLIFGLGAVWGLGVDPSLQIDVLTLEL  
GESTTGVRQVPLHNGTKAFLFQDTPRSIKASTATAE  
QFFQKLRNKHEFTILVTLKQTHLNSGVLSIHHLDRYL  
ELESSGHRNEVRLHYRSGSHRPHEVFPYILADDKWH  
KLSLAISASHLILHIDCNKIYERVVEKPSTDLPLGTTFWL  
GQRNNAHG YFKGIMQDVQLLVMPOGFIAQCPDLNR  
TCPTCNDFHGLVQKIMELQDILAKTSAKLSRAEQRMN  
RLDQCYCERTCTMKGTTYREFESWIDGCKNCTCLNG  
TIQCETLICPNPDCPLKSALAYVDGKCKECKSICQF  
QGRTYFEGERNVYSSSGVCLYECKDQTMKLVES  
GCPALDCPESHQITLSHSCCKVCKGYDFCSERHNCM  
ENSICRNLDRAVCSRDGFALREDNAYCEDIDECA  
EGRHYCRENTMVCVNTPGSFMCICTGYIRIDDYSCTE  
HDECITNQHNCDENALCFNTVGGHNCVCKPGYTG  
GTTCKAFCKDGCRRNGGACIAANVCACPGFTGPSC  
ETDIDECSDGFVQCDSRANCINLPGWYHCECRDGYH  
DNGMFSPSGESCEDIDECGTGRHSCANDTICFNLDG  
GYDCRCPHGKNCIGDCIHDGKVKHNGQIWWLEND  
CSVCSQNGFVMCRRMVDCENPTVDFCCPECDP  
RLSSQCLHQNGETLYNSGDTWWQNCQQCRCLQGEV  
DCWPLPCPDVECEFSILPENECCPRCVTDPCQADTIR  
NDITKTCLEDMNVVRF TGSSWIKHGTECLCQCKNG  
HICCSVDPQCLQEL  
(1 - 816 aa encoded by BC020544)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W: ptglab.com  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.