

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

XIAP Fusion Protein



Basic Information

Catalog Number:

Ag0725

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MTFNSFEQSKTQVPADINKEEEFVEEFNRLKTFANFPS
GSPVSASTLARAGFLYTGEDTVRCFSCHAAVDRWQ
YGDSAVGRHRKVPNCRFINGFYLENSATQSTNSGIQ
NGQYKVENYLGRDHFALDRPSETHADYLLRTGQVV
DISDTIYPRNPAMYSEEARLKSFNWPDYAHLPRELA
SAGLYYTGIGDQVQCFCCGGKLNWEPDRAWSEH
RRHFPNCFVVLGRNLNIRSEDAVSSDRNFPNSTNLPR
NPSMADYEARIFFGTWIVSNKEQLARAGFYALGEG
DKVKCFHCGGGLTDWKPSEDPWEQHAQWYPGCKYL
LEQKQGEYINNIHLTHSLEECVLRTEKTPSLTRRIDDI
FQNPVMQEAIRMGFSFKDIKKIMEEKIQISGSNYKSLE
VLVADLVNAQKDSMQDESSQTSLQKEISTEEQLRRL
QEEKLCIKMDRNIAIVFVPCGHLVTCCKQCAEAVDK
CPMCYTVITFKQKIFMS

(1-497 aa encoded by BC032729)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.