

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant CYP302a1 protein



Basic Information

Catalog Number:

Ag0761

Size:

50 µg

Form:

Available lyophilized

Species:

Expression Source:

e coli-derived, PET28a, with N-terminal 6^{His}.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MYLSPVKSIKSATLMTSRSCATVVLENVKPNQIPGP
RGPFGLGNLYQYIPGIGKYSFDALHESGQDKYEKYGPI
VRETMVPGQDIVWLYDPNDIAAVLNDKTPGIYPSRRS
HTALAKYRRDRPNVYRTAGLLATNGIEWWKIRSELQK
GLSSPQSVRNFLPLTDKVTREFVASMNSTEHDCVPDF
MPAISRLNLELICVMAFDVRLDSFSDEQMKPNSLSSRL
MESAEVTNQSIPLTDQGFQLWKYFETPAYRKLKKAQE
FMEKTAVELVSQKLLYFDEDQQLASGRHRSRLLEEY
LRNPNLELHDIIGMAADLLLAGVHTSSYTAFALYHLCL
NPDAQDKLYQEACRILPDPWECQIEAAALNSEASYC
RAVLKESLRLNPISIGVGRILNKDATLGGYHVPKGTVV
VTQNLVSCRQERYFKNPTKFI PERWMRETKEVDNPNYL
VLPGHGMRSCLARRMAEQNMLVLLLRILRSYEIDWKG
KVPMNIETKLINQPDQPIKIAFRSRKS
(1-515 aa encoded by AY947549)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.