

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# WNT7A Fusion Protein



## Basic Information

**Catalog Number:**  
Ag0874

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Validated Application:**  
Blocking peptide

**Peptide Sequence:**  
MNRKARRCLGHLFSLGMVYLRIIGFSSVVALGASIIIC  
NKIPGLAPRQRAICQSRPDIIIVIGESQMGLDECQF  
QFRNGRWNCALGERTVFGKELKVGSRFAAFTYAIIA  
AGVAHAITAACTQGNLSDCGCDKEKQGOYHRDEG  
WKWGGCSADIRYIGFAKVFVDAREIKQNARTLMNL  
HNNEAGRKILEENMKLECKCHGVSGSCTTKCWTTLP  
QFRELGYVLKDKYNEAVHVEPVRASRNKRPTFLKIKKP  
LSYRKPMDDLVIYEKSPNYCEEDPVTGSGVTQGRAC  
NKTAPQASGCDLMCCGRGYNTHQYARVWQCNCCKF  
HWCCYVKCNTCSERTEMYTCK  
(1-349 aa encoded by BC008811)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) E: proteintech@ptglab.com W:  
(toll free in USA), or 1(312) 455-8498 ptglab.com  
(outside USA)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.