

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human CRLF1 protein



Basic Information

Catalog Number:

Ag10396

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

HTAVISQDPTLLIGSSLLATCSVHGDPGATAEGLY
WTLNGRRLPPELSRVLNASTLALANLNGSRQ RSGD
NLVCHARDGSILAGSCLYVGLPPEKPVNISCWSKNMK
DLTCRWTPGAHGETFLHTNYSKYLKRWYGDNTCE
EYHTVGPHSCHIPKDLALFPTYEIWWEATNRLGSARSD
VLTLDILDVVTTDPPPDVHVS RVGGLEDQLSVRWVSP
PALKDFLFQAKYQIRYRVEDSVDWKVVDDVSNQTSC
RLAGLKPGTVYFVQVRCNPFGIYGSKKAGIWEWSHP
TAASTPRSERPGPGGGACEPRGGEPSSG PVRRELKQF
LGWLKKAHAYCSNLSFRLYDQWRAWMQKSHKTRNQD
EGILPSGRRGTARGPAR

(39-422 aa encoded by BC044634)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.