

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# ACAD10 Fusion Protein



## Basic Information

**Catalog Number:**

Ag10529

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MGGVLI P S P G R V A A E W E V Q N R I P S G T I L K A L M E G G E N  
G P W M R F M R A E I T A E G F L R E F G R L C S E M L K T S V P V D S F  
F S L L T S E R V A K Q F P V M T E A I T Q I R A K G L Q T A V L S N N F Y  
L P N Q K S F L P L D R K Q F D V I V E S C M E G I C K P D P R I Y K L C L E  
Q L G L Q P S E S I F L D D L G T N L K E A A R L G I H T I K V N D P E T A  
V K E L E A L L G F T L R V G V P N T R P V K K T M E I P K D S L Q K Y L K  
D L L G I Q T T G L

(1-237 aa encoded by BC015056)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: proteintech@ptglab.com  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.