

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human ZNF3 protein



## Basic Information

### Catalog Number:

Ag10578

### Size:

50 µg

### Form:

Available lyophilized

### Species:

human

### Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

### Biological Activity:

Not tested

### Endotoxin Level:

Please contact the lab for more information

### Peptide Sequence:

METQADLVSQEPQALLDSALPSKVPAFSDKDSLGD  
MLAAALLKAKSQELVTFEDVAVYFIRKEWKRLPAQR  
DLYRDVMLENYGNVFSLDRETRTENDQEISEDTRSHG  
VLLGRFQKDISQGLKFEAYEREVSLKRPLGNSPGERL  
NRKMPDFGQVTVEEKLTPRGERSEKYNDFGNSFTVNS  
NLISHQRLPVGDRPHKCDCECKSFNRTSDLIQHQRHT  
GEKPYECNECGKAFSQSSHLIQHQRHTGEKPYECSD  
CGKTFSCSSALILH

(1-274 aa encoded by BC013603)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.