

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human ULK4 protein



## Basic Information

Catalog Number:

Ag10784

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

I L C E D P L P P I P K D S S R P K A S S D F I N L L D G L L Q R D P Q K R L T  
W T R L L Q H S F W K K A F A G A D Q E S S V E D L S R N T M E C S  
G P Q D S K E L L Q N S Q S R Q A K G H K S G Q P L G H S F R L E N P T  
E F R P K S T L E G Q L N E S M F L L S S R P T P R T S T A V E V S P G E D  
M T H C S P Q K T S P L T K I T S G H L S Q Q D L E S Q M R E L I Y T D S  
D L V V T P I I D N P K I M K Q P P V K F D A K I L H L P T Y S V D K L L F L  
K D Q D W N D F L Q Q V C S Q I D S T E K S M G A S R A K L N L L C Y L  
C V V A G H Q E V A T R L L H S P L F Q L L I Q H L R I A P N W D I R A K  
V A H V I G L L A S H T A E L Q E N T P V V E T T S S I G I L N C L V Q  
H S T P V P R Q C L V Y V

(231-580aa encoded by BC014794)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.