

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human CHDH protein



Basic Information

Catalog Number:

Ag10883

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

Y L H P A L S R T N L K A E A E T L V S R V L F E G T R A V G V E Y V K N G
Q S H R A Y A S K E V I L S G G A I N S P Q L L M L S G I N A D D L K K L
G I P V V C H L P G V G Q N L Q D H L E I Y I Q Q A C T R P I T L H S A Q
K P L R K V C I G L E W L W K F T G E G A T A H L E T G G F I R S Q P G V
P H P D I Q F H F L P S Q V I D H G R V P T Q Q E A Y Q V H V G P M R
G T S V G W L K L R S A N P Q D H P V I Q P N Y L S T E T D I E D F R L C
V K L T R E I F A Q E A L A P F R G K E L Q P G S H I Q S D K E I D A F V R
A K A D S A Y H P S C T C K M G Q P S D P T A V V D P Q T R V L G V E
N L R V D A S I M P S M V S G N L N A P T I M I A E K A A D I I K G Q P
A L W D K D V P V Y K P R T L A T Q R

(244-594aa encoded by BC034502)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.