

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

HLA-DQB2 Fusion Protein



Basic Information

Catalog Number:

Ag11284

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MSWKMALQIPGGFWAAAVTVMLVMLSTPVAEARDF
PKDFLVQFKGMCYFTNGTERVIRGVARYIYNREEYGRF
DSDVGEFQAVTELGRSIEDWNNYKDFLEQERAAVDK
VCRHNYEAE LR T T L Q R Q V E P T V T I S P R T E A L N H H N L L
V C S V T D F Y P A Q I K V Q W F R N D Q E E T A G V V S T S L I R N G
D W T F Q I L V M L E I T P Q R G D I Y T C Q V E H P S L Q S P I T V E W
R P R G P P P A G L L H

(1-231 aa encoded by BC031995)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.