

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human GPBP1L1 protein



## Basic Information

**Catalog Number:**

Ag11316

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

MAQHDFVPAWLNFTPQSAKSPTATFEKHGEHLPRG  
EGRFGVSRRRHNSDGFNFNGPLRTAGDSWHQPSLF  
RHDSVD SGVSKGAYAGITGNPSGWSSSRGHDGMS  
QRSGGGTGNHRHWNGSFHSRKGCAFQEKPPMEIREE  
KKEDKVEKLFEEEDFPSLNPEAGKHQPCRPIGTPS  
GVWENPPSAKQPSKMLVIKVKSKEDPAAAFSAAFS  
PGSHHANGNKLSSVVPVYKLVKPVPPSPKPNAW  
KANRMEHKSGSLSSRESAFTSPISVTKPVVLSGAAL  
SSPKESPSSTPPIEISSRLTKLRRTRDRKSEFLKTLKD  
DRNGDFSENRC DKLEDLEDNSTPEPKENGEEGCHQ  
NGLALPVVEEGEVLSHSLEAEHRLKAMGWQYEPEND  
ENCLPLTEDELKEFHMKTEQLRRNGFGKNGFLQSRSS  
LFSWRSTCKAEFEDSDTETSSETSDDDAWK  
(1-474 aa encoded by BC022218)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.