

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human MKX protein



Basic Information

Catalog Number:

Ag11362

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MNTIVFNKLSSAVLFEDGGASERERGGRPYSGVLDSP
HARPEVGIPDGPPLKDNLGLRHRRTGARQNGGKVRHK
RQALQDMARPLKQWLYKHRDNPYPTKEKILLALGSQ
MTLVQVSNWFANARRRLKNTVRQPDLSWALRIKLYN
KYVQGNAERLSVSSDDSCSEDGENPPRTHMNEGGYN
TPVHHPIKSENSVIKAGVRPESRASEDYVAPPKYKSS
LLNRYLNDSLRHVMA T N T T M M G K T R Q R N H S G S F S S N
E F E E E L V S P S S E T E G N F V Y R T D T L E N G S N K G E S A A N R
K G P S K D D T Y W K E I N A A M A L T N L A Q G K D K L Q G T T S C I I
Q K S S H I A E V K T V K V P L V Q Q F

(1-352aa encoded by BC036207)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)
(toll free in USA), or 1(312) 455-8498
(outside USA)

E : proteintech@ptglab.com W:
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.