

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human ASTN2 protein



## Basic Information

Catalog Number:

Ag11409

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MNTLLCKGMFCLLSWEADSRGRLGEYTLQPLSLQTEE  
TTELGSKKELKSMFFITYLSGLLTAQMLSDDQLISGVEI  
RCEEKGRCPSTCHLCRRPGKEQLSPTVLEINRVVPL  
YTLIQDNGTKEAFKALMSSYWCSGKGDVIDDWCRC  
DLSAFDANGLPNCSPLLQPVLRLSPTVEPSSTVVSLEW  
VDVQPAIGTKVSDYILQHKKVDEYTDLTLYGFLSF  
ADDLLSGLGTSCVAAGRSHGEVPEVSIYVIFKCLEPD  
GLYKFTLYAVDTRGRHSELSTVTLRTACPLVDDNKAAE  
IADKIYNLYNGYTSGKEQQMAYNTLMEVSASMLFRV  
QHHYNSHYEKF GDFVWRSEDELGPRKAHLIRRLERV  
SSHCSSLLRSAYIQSRVETVPYLF CRSEEV RPAGMVV  
YSILKDTKIMCEEKMMVSMARNTYGESKGR  
(1-440 aa encoded by BC029272)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.